


## TRAYS




- Precision, pressure-formed PETG
- 1073 B Tyvek Lidstock
- Single or double barrier
- Custom foam and urethane protection available
- Five year shelf life (Data Upon Request)
- EO and Gamma sterilization compatible

\*NOTE: All dimensions are approximate

Tray #	Inside Dimensions			Outside Dimensions		
	L	W	D	L	W	D
RP1MD	4.5"	4.5"	2.8"	5.5"	5.5"	2.9"
RP1LG	5.6"	5.6"	3.2"	6.4"	6.4"	3.3"
RP2SM	3.4"	3.4"	2.4"	4.6"	4.6"	2.5"
RP2MD	4.5"	4.5"	2.8"	5.5"	5.5"	2.9"
RP2.5SM	3.4"	3.4"	1.2"	4.6"	4.6"	1.3"
RP2.5MD	4.5"	4.5"	1.4"	5.5"	5.5"	1.5"
RP2.5XS	2.3"	2.3"	1.0"	3.7"	3.7"	1.1"
RP3MD	9.8"	4.1"	1.2"	11.1"	5.2"	1.3"
RP3LG	10.0"	4.3"	1.4"	12.3"	6.4"	1.5"

## POUCHES



- Chevron, peelable structure
- Materials: Tyvek 1073B/polyester, nylon, or foil
- Single or double barrier
- Five year shelf-life
- EO (Tyvek only) and gamma compatible

\*NOTE: All dimensions are approximate


Pouch #	Inside Dimensions		Outside Dimensions	
	L	W	L	W
1097	8.5"	1.75"	9.25"	2.5"
1098	8.5"	2.75"	9.25"	3.5"
1099	2.0"	1.25"	4.25"	2.0"
1100	4.0"	2.25"	7.0"	3.0"
1101	5.0"	3.25"	9.0"	4.0"
1102	6.0"	5.25"	10.0"	6.0"
1103	14.0"	11.25"	18.0"	12.0"
1104	6.25"	2.25"	7.0"	3.0"
1105	4.25"	1.25"	5.0"	2.0"
1106	4.5"	2.75"	5.25"	3.5"
1107	6.25"	3.75"	7.0"	4.5"

## FOAM/URETHANE

- Die-cut and molded
- Cross-linked polyethylene foam and urethane
- EO and gamma sterilization compatible
- Custom designs

\*NOTE: All dimensions are approximate

#	Dimensions		
	L	W	D
A-Foam	4.4"	4.4"	2.0"
B-Foam	3.4"	3.4"	0.13"
C-Foam	4.8"	4.8"	0.38"
C-Foam	2.3"	2.3"	0.13"
E-Foam	6.5"	6.5"	0.13"
F-Foam	1.5"	2.0"	N/A
G-Uret	9.8"	3.9"	N/A
H-Foam	16.50"	8.0"	2.0"
I-Foam	1.13"	1.13"	N/A



## CARTONS

- SBS and e-flute corrugated materials
- Printed and unprinted
- Custom graphic designs
- Custom sizes

\*NOTE: All dimensions are approximate

#	Dimensions		
	L	W	D
1001	6.8"	3.2"	6.8"
1002	4.375"	1.312"	6.0"
1003	6.0"	1.3"	6.0"
1004	2.8"	0.8"	7.0"
1005	5.8"	3.0"	5.8"
1006	6.5"	1.6"	12.5"
1007	3.0"	0.875"	2.50"
1008	9.0"	4.0"	7.0"
1010	6.0"	1.0"	13.0"

